



PMI Feeds, Inc.®

PMI® Micro-Stabilized Rodent Liquid Diet LD 101/101A*

Description Diet LD 101A is a dry powder used to prepare a liquid diet for rodents in alcohol studies. The powder is designed to be mixed with alcohol and carbohydrate prior to feeding. When mixed according to instructions, it provides a similar level of nutrition as LD 101.

Diet LD 101 is a dry powder used to prepare a liquid diet formulated for rodents. It is nutritionally balanced with excellent palatability. This diet is especially designed for use as the control diet when LD 101A is being used as the test diet. LD 101 can be used in other applications where solid diets are not appropriate.

Various flavors are available upon request.

Features and Benefits

- Nutritionally-balanced
- Volatile ingredients can be included
- Easily prepared
- Provides stable nutrients
- Shipped in dry form to simplify storage, shipping and stability
- Minimal foaming
- Fully suspended
- Stabilized against microbial growth

Typical Analysis

Crude protein not less than.....	3.8%
Crude fat not less than.....	3.9%
Crude fiber not more than.....	1.0%
Ash not more than.....	0.4%

Values are based upon the liquid form of the diet when prepared according to directions.

***Diet Preparation Instructions:** To 770 gms. deionized/distilled water, add 230 gms. Micro-Stabilized Rodent Liquid Diet mix (LD 101). Blend vigorously for 15-30 seconds with a mechanical blender until completely suspended. For best results add water to blender before dry mix. See information provided to prepare isocaloric diets with alcohol or other test substances.

Ingredients

Vitamin-free casein, olive oil, maltodextrin, dried corn syrup, soy fiber, corn oil, suspension colloid, safflower oil, L-cystine, DL-methionine, vitamin A acetate, cholecalciferol, DL-alpha tocopheryl acetate, menadione sodium bisulfite (source of vitamin K), fumaric acid, citric acid, propionic acid, ascorbic acid, potassium sorbate, cyanocobalamin, thiamin mononitrate, riboflavin, calcium pantothenate, nicotinic

acid, choline chloride, pyridoxine hydrochloride, folic acid, inositol, p-aminobenzoic acid, biotin, calcium acetate, calcium phosphate, potassium phosphate, sodium phosphate, magnesium sulfate, sodium chloride, manganese sulfate, ferrous fumarate, zinc chloride, cupric sulfate, chromium chloride, sodium fluoride, ammonium molybdate, calcium iodate and sodium selenite.

Feeding Directions

Diet consumption will vary according to animal size and sex. An average rat should consume at least 50 grams of diet daily to sustain a creditable average daily weight gain. The growth rate of rats maintained on this diet should be similar to that attained by young rats (55-100 gram) maintained on a good quality, nonpurified rodent diet. Mice should consume at least 20 grams per day.

Allow new animals an adequate period of time to adjust to their surroundings. After they have adjusted, introduce the liquid diet gradually by offering some of the liquid diet while the regular diet is still present. Gradually decrease the amount of regular diet offered

while increasing the amount of liquid diet over a 3-5 day period. Additional time for adjustment may be necessary for the ethanol diets.

Prepare the diet as frequently as needed and always refrigerate to minimize loss of nutrients. Fresh diet should be prepared at least every 5 days. Although the diet may be bacteriologically sound for a longer period of time, diet more than 5 days old may have deteriorated nutritionally. Before using diet which has been prepared on a previous day, check to ensure all of the ingredients are in suspension. Remix if necessary. Additional water may be provided in separate drinking tubes, but may not be consumed.



PMI Feeds, Inc.®

PMI® Micro-Stabilized Rodent Liquid Diet LD 101/101A*

Diet Preparation Instructions

PMI® Micro-Stabilized Rodent Liquid Diet (LD 101)

Water (deionized/distilled)..... 770 grams
Micro-Stabilized Rodent Liquid Diet
mix (LD 101)..... 230 grams

Blend vigorously for 15-30 seconds with a mechanical blender until completely suspended. For best results add water to blender first, then add dry mix.

PMI® Micro-Stabilized Alcohol Rodent Liquid Diet (LD 101A)

Diet composition varies according to the amount of alcohol added to maintain an isocaloric diet. The following chart indicates the amounts of deionized/distilled water, PMI® Micro-Stabilized Alcohol Rodent Liquid Diet LD 101A mix (Dry Mix), PMI® Maltodextrin LD 104, and ethanol to be used to make one kilogram of liquid diet.

<u>% Calories from Ethanol</u>	<u>gms. Water</u>	<u>gms. Dry Mix</u>	<u>gms. Maltodextrin</u>	<u>gms. Ethanol</u>
36	810	140	0	50
30	803.7	140	14	42.3
20	792.8	140	39	28.2
10	783	140	63	14.1

Add water to a mechanical blender, then add dry mix and maltodextrin. Blend vigorously for 15 seconds. Add ethanol to the blended mixture and blend for an additional 15 seconds.

For calculation purposes:

- 140 gms. dry Alcohol Rodent Liquid Diet mix = 645 kcal.
- Ethanol = 7.1 kcal/gm
- PMI® Maltodextrin LD 104 = 3.96 kcal/gm

Additional Considerations:

- For best results a mechanical blender should be used for diet preparation; hand blending does not suspend the diet adequately to avoid some settling out of undissolved ingredients.
- Do not over-blend; excessive mechanical blending creates foaming.
- Prepare fresh diet at least every **5 days**. Keep unused portion of diet refrigerated or at room temperature (not to exceed 75°F). Before using diet which has been prepared on a previous day, check the diet to insure all of the ingredients are still in suspension - remix if necessary. Although the diet will remain bacteriologically stable for up to 10 days, product more than **5 days** old should not be used due to changes in the nutritional quality.

* Product Code Numbers



PMI Feeds, Inc.

LabDiet
THE RICHMOND STANDARD

PMI® Micro-Stabilized Rodent Liquid Diet LD 101/101A*

Chemical Composition¹

Nutrients²

Protein gm/kg	38.4
Arginine gm/kg	1.38
Cystine gm/kg	0.24
Glycine gm/kg	0.76
Histidine gm/kg	1.01
Isoleucine gm/kg	1.88
Leucine gm/kg	3.4
Lysine gm/kg	2.86
Methionine gm/kg	1.21
Phenylalanine gm/kg	1.88
Tyrosine gm/kg	1.99
Threonine gm/kg	1.52
Tryptophan gm/kg	0.43
Valine gm/kg	2.24
Fat gm/kg	39.6
Fiber (Crude) , gm/kg	10.0

Minerals

Calcium mg/kg	1300
Phosphorus (total) mg/kg	1000
Potassium mg/kg	875
Magnesium mg/kg	125
Sulfate mg/kg	500
Sodium mg/kg	323
Chlorine mg/kg	390
Fluorine, mg/kg	0.25
Iron, mg/kg	17.6
Zinc, mg/kg	7.6
Manganese, mg/kg	13.7
Copper, mg/kg	1.5
Chromium, mg/kg	0.6
Iodine, mg/kg	0.05
Molybdenum, mg/kg	0.11
Selenium, mg/kg	0.027

Vitamins

Menadione, mg/kg	0.24
Thiamin mononitrate, mg/kg	1.6
Riboflavin, mg/kg	1.5
Nicotinic acid, mg/kg	7.5
Calcium pantothenate, mg/kg	4.0
Choline, mg/kg	250
Folic Acid, mg/kg	0.5
Pyridoxine hydrochloride, mg/kg	1.75
Biotin, mg/kg	0.05
Inositol, mg/kg	25.0
p-aminobenzoic acid, mg/kg	12.5
B ₁₂ , mg/kg	0.025
Vitamin A, IU/kg	3000
Vitamin D ₃ (added), IU/kg	400
Vitamin E, IU/kg	30
Ascorbic Acid, mg/gm	10.0

*Product Code Numbers

¹ Values are based upon the liquid form of the diet when prepared according to directions. Based on the latest ingredient analysis information.

Energy*

Protein, kcal/kg	164
Fat, kcal/kg	350
Carbohydrates, kcal/kg	486

*Energy Levels used (kcal/gm)

Protein = 4.27; Fat = 8.84; Carbohydrates = 3.96; Ethanol = 7.1. These values are different than the 4, 9, 4 kcal/gm for protein, fat & carbohydrate, respectively, as generally used.

* Lieber, CS & LM DeCarli (1982) Alcoholism: Clinical and Experimental Research 6: 523-531.
Miller, SS, ME Goldman, CK Erickson & RL Shorey (1980) Psychopharmacology 68: 55-59.

* 1 kilogram of diet in liquid form, when prepared according to directions, provides 1000 kilocalories (1 kcal per gram).

KODAK MEGAPLUS Camera,

Models:

ES 1.0

ES 1.0/10 Bit

ES 1.0/TH

- I** Optomechanical Specification
- II** Imaging Performance Specification

REVISION 9.0

KODAK and MEGAPLUS are trademarks.
All other trademarked names are the property of their respective companies.

TABLE OF CONTENTS

<u>Section</u>	<u>Sub-section</u>
I	OPTOMECHANICAL
	Image sensor 1
	Sensor cover glass 2
	Optical window 3
	Electronic Focal plane shutter 4
	Lens interface 5
	Drawings 6
II	IMAGING PERFORMANCE
	Signal Output 1
	Responsivity and Equivalent ISO speed 2
	Gamma ratio 3
	Dynamic range 4
	Maximum Luminous Signal/Noise ratio 5
	Sensor defects and classes 6
	Dark field noise and nonuniformity 7
	Nonlinearity 8
	Minimum irradiance and illuminance 9
	Photoresponse nonuniformity 10
	PHOTON TRANSFER FUNCTION and GAMMA TRANSFER FUNCTION 11
	MODULATION TRANSFER FUNCTION 12
	Limiting Resolution 13
	CONTRAST TRANSFER FUNCTION 14
	Spectral responsivity 15
	Crosstalk and Smear 16
	Blooming Suppression 17
	Stray light DATA NOT AVAILABLE 18
	Valid signal output and measurement accuracy 19
	Measurement instrumentation 20

I OPTO-MECHANICAL SPECIFICATION

1) Image sensor

MODEL	Kodak KAI 1001M CCD		
TYPE	Front illuminated interline architecture		
PHOTOSENSITIVE AREA	9.1 H by 9.2 V	mm's	
ASPECT RATIO	1:1	ratio	
HORIZONTAL PIXELS	1008	valid pixels	
VERTICAL PIXELS	1018	valid pixels	
PIXEL DIMENSIONS	9 H by 9 V	microns	
FILL FACTOR	60	percent	

DEFECTS refer to PART II, SECTION 7

2) Sensor cover glass

GLASS TYPE	Corning 7059		
REFRACTIVE INDEX	1.5	@	587.6 nm
THICKNESS	0.76	mm +/-	0.05 mm
SURFACE FINISH	5	microns	
COATING	Anti-reflective	(both sides)	optimized 430-630 nm's
REAR SURFACE LOCATION	2.5	mm	in front of the sensor plane

3) Optical window

CLEAR WINDOW

GLASS TYPE	BK7		
REFRACTIVE INDEX	1.52	@	587.6 nm
DIAMETER	31.7	mm +/-	0.3 mm
THICKNESS	1.00	mm +/-	0.05 mm
FLATNESS	0.2	waves @	546.1 nm over 80% of surface
SURFACE FINISH	60-40	Scratch and Dig	
OPTIONAL COATINGS	Wide band hot mirror		

INFRARED FILTER

GLASS TYPE	Schott KG5 filter glass		
REFRACTIVE INDEX	1.51	@	587.6 nm
DIAMETER	31.7	mm +/-	0.2 mm
THICKNESS	1.00	mm +/-	0.1 mm
FLATNESS	2-3	Waves @	546.1 nm
SURFACE FINISH	60-40	Scratch and Dig	
COATING	None		

4) Focal plane shutter

TYPE	fully electronic	
RESOLUTION	0.0000625	seconds
CLEAR APERTURE	9.1 H by 9.2 V	mm's
EXPOSURE RANGE	0.000125 and up	seconds

The shutter may exhibit less than 100% repeatability from frame-to-frame, as well as some spatial nonuniformity of exposure within any single frame. Repeatability and nonuniformity are reported as percent of the mean gray level output.

DEFINITION:	Nonuniformity Defined as the spatial nonuniformity of exposure within any single frame.
DEFINITION:	Repeatability Defined as the variation in the mean signal output level for any single frame compared to any other.
DEFINITION:	Lifetime Defined in terms of either the mean number of continuous shutter operating hours or number of shutter exposure cycles before complete mechanical failure of the shutter mechanism.

NONUNIFORMITY	100	percent	@ <= 0.006 seconds exposure
REPEATABILITY	NA	percent	
LIFETIME	NA	cycles	

NOTE:	The rated <i>lifetime</i> and <i>repeatability</i> of camera focal plane shutters applies to mechanical shutters only. <i>Lifetime</i> is based on complete mechanical failure and does not consider degradation of performance with age. Nonuniformity and repeatability may degrade with age.
	Repeatability and nonuniformity of mechanical shutters may degrade with exposure times < 0.050 seconds.

Measurement conditions:

- Infrared-filtered Xenon light source*
- Integrating sphere illuminator*
- 25 ° C ambient operating temperature*
- Digital video output*
- Defect concealment disabled*
- Center and edge of frame ROI image samples*

5) Lens interface

TYPE	C mount		
FLANGE FOCAL DISTANCE	17.5	mm	+/- 0.05
NONCOPLANARITY	0	arc-minutes	+/- 47.4
DECENTER	0	mm	+/- 0.45

6) Drawings





FIGURE 1.1


ES 1.0 and ES 1.0/10 BIT



SIDE VIEW

FIGURE 1.2


ES 1.0 and ES 1.0/10 BIT



BOTTOM VIEW

FIGURE 1.3


ES 1.0/TH



SIDE VIEW

FIGURE 2.3


ES 1.0/TH



SIDE VIEW

FIGURE 2.4

ES 1.0/TH



FRONT VIEW

FIGURE 3


II IMAGING PERFORMANCE SPECIFICATION

The performance data contained in this document is presented as "typical" and does not represent a guarantee of any specific performance level. Where noted, some performance data is reported as theoretical.

1) Signal Output

TYPE	Digital RS 422 Differential		
FORMAT	8 or 10 Bit parallel		
FRAME RATE	MINIMUM	MAXIMUM	
Continuous	30	30	frames/second
Triggered	0	30	frames/second
Controlled	0	30	frames/second
Triggered Double Exposure	0	15	frames/second

The ES 1.0 internally digitizes video using 10 significant bits of gray scale resolution. The Camera will then output either a the full 10 bit data or a select 8 bit set of the 10 bit internal representation. The ES 1.0 *base model* output is limited to the select eight ('8') significant bits of digital video. Invoking the Camera's 'DIGITAL GAIN' command selects the specific 8 bits for output as follows:



		<u>8 BIT Output</u>	<u>10 BIT Output</u>
-	Command 'DGN 1'	Bits 2-9	Bits 0-9
-	Command 'DGN 2'	Bits 1-8	
-	Command 'DGN 4'	Bits 0-7	

2) Responsivity and Equivalent ISO speed

DEFINITION: Responsivity
 This is defined as the camera's nominal response to light. In other words, the change in the camera's digital video output level, measured in counts, for a given change in illuminance or irradiance delivered to the camera. Radiometric responsivity is the derivative of the GAMMA TRANSFER FUNCTION measured @ 550 nm's. Photopic responsivity is the derivative of the GAMMA TRANSFER FUNCTION measured using a specified broad band illumination.

DEFINITION: Equivalent ISO speed
 Equivalent ISO speed is a variation on the ISO photographic speed standard used to specify the 'sensitivity' of many film-based imaging systems. $\text{Equivalent ISO speed} = \frac{78}{H_{Sat}}$ where H_{Sat} = focal plane exposure in Lux-seconds for the maximum VALID digital video output count. Refer to ISO STANDARD #12232, "Photography - Electronic Still Picture Cameras - Determination of ISO Speed"

3) Gamma ratio

DEFINITION: Gamma ratio
 Defined as the derivative of the best fit line to a 'Log transformed' GAMMA TRANSFER FUNCTION. The GAMMA TRANSFER FUNCTION is 'Log transformed' by applying a $\text{Log}_{\text{Base } 10}$ remapping function to both the X and Y coordinate values of the measured GAMMA TRANSFER FUNCTION. See GAMMA TRANSFER FUNCTION.
 Linear GAMMA TRANSFER FUNCTION \Rightarrow $\text{OUTPUT} = m * (\text{INPUT}) + \text{offset}$
 $m = \text{responsivity,}$
 $\text{offset} = \text{black background level}$
 Gamma ratio \Rightarrow $\text{Log}[\text{OUTPUT}]/\text{Log}[m * \text{INPUT}] + \text{offset}$
 Log Transformed GAMMA TRANSFER FUNCTION \Rightarrow $\text{Log}[\text{OUTPUT}] = \text{Gamma ratio} * \text{Log}[m * \text{INPUT}] + \text{offset}$
 or simply \Rightarrow $\text{OUTPUT} = m * \text{INPUT}^{\text{Gamma ratio}} + \text{offset}$

<u>GAIN</u> (DGN 'X')	<u>Gamma ratio</u> (8 bit)	<u>Gamma ratio</u> (10 bit)	
1	1	1	ratio
2	1	-	"
4	1	-	"

4) Dynamic Range

DEFINITION: Dynamic Range
 Defined on the PHOTON TRANSFER FUNCTION as the ratio between the $\text{SQRT}[\text{Y axis coordinate value}]$ (i.e. signal variance axis) for the lowest VALID mean signal output level, to the X axis coordinate value (i.e. signal variance axis) for highest VALID mean signal output; expressed in Decibels as $20 * \text{LOG}_{10}(\text{Ratio})$.

<u>GAIN</u> (DGN 'X')	<u>Dynamic Range</u> (8 bit)	<u>Dynamic Range</u> (10 bit)	
1	56	59	Decibels
2	52	NA	"
4	46	NA	"

5) Maximum Luminous Signal/noise ratio

DEFINITION: Maximum Luminous Signal/noise ratio
 Defined on the PHOTON TRANSFER FUNCTION as the ratio between the $\text{SQRT}[\text{Y axis coordinate value}]$ (i.e. signal variance axis) for the highest VALID mean signal output level (i.e. near A/D saturation), to the X axis coordinate value (i.e. mean signal level) for highest VALID mean signal output level (i.e. near A/D saturation); expressed in Decibels as $20 * \text{LOG}_{10}(\text{Ratio})$.

<u>GAIN</u> (DGN 'X')	<u>Maximum Luminous Signal/Noise Ratio</u> (8 bit)	<u>Maximum Luminous Signal/Noise Ratio</u> (10 bit)	
1	46	46	Decibels
2	42	NA	"
4	39	NA	"

6) Sensor defects and classes

The Camera can be equipped with an image sensor of any of three general quality levels. General sensor quality is specified with regard to the number and location of pixels which exhibit *markedly different output* from their neighbors when operated under conditions specified below. These nonuniformities are called *defects*. Defects can be of the *point, cluster, or column/row* type.

DEFINITION:	Saturation
This is defined as the pixel illumination level wherein no further increase in camera output is possible and either vertical smearing or blooming suppression begin.	
DEFINITION:	Point defect
An isolated defective pixel.	
DEFINITION:	Major Defective Pixel
This is defined as a pixel whose signal deviates by more than 7.5% from the mean value of all active pixels in complete darkness, or by more than 15% from the mean value of all active pixels under uniform illumination at 80% of saturation.	
DEFINITION:	Minor Defective Pixel
This is defined as a pixel whose signal deviates by more than 8 millivolts from the mean value of all pixels in complete darkness.	
DEFINITION:	Cluster defect
A grouping of 2 to 4 contiguous major defective pixels.	
DEFINITION:	Column or row defect
A grouping of more than 4 contiguous major defective pixels along a single column or row.	

CLASS	Point Defects		Cluster Defects	Column/Row Defects
	MAJOR	MINOR		
CLASS 0	<= 0	<= 50	0	0
CLASS 1	<= 5	<= 50	0	0
CLASS 2	<= 20	<= 100	<= 4	<= 0
CLASS 3	<= 50	<= 400	<= 8	<= 4

OPERATING CONDITIONS:		
Illumination level	80	% of Saturation
Sensor operating Temperature	40	Degrees Celsius
Exposure time	0.070	seconds
Readout time	0.070	seconds
Gain	0	Decibels

For the Model ES 1.0 cameras, 80% of saturation is reached when the Camera is illuminated to approximately full scale output (full scale = 255 or 1023 for 8 or 10 bits respectively) for DGN 1 gain setting.

NOTE: Nonuniformities which represent known sensor defects have been pre-identified using only the above specified sensor operating conditions. The Model ES 1.0 cameras have a typical sensor operating temperature of 40° Celsius. Higher sensor operating temperature may reveal additional nonuniformities not called out in the defect map which accompanies the Camera. Higher camera gain settings, longer exposure times, and/or higher or lower illumination levels may also in some cases increase the number of dark field nonuniformities or photoresponse nonuniformities that will qualify as defects according to the above defect definitions.

7) Dark field noise and nonuniformities

DEFINITION: Pixel-to-pixel dark field noise
This is a random pixel-to-pixel video output level disparity within a single frame and from frame-to-frame. It is sourced in the Camera's analog channel electronics and individual pixel dark current variations. It increases with increasing gain, exposure time, and increasing ambient operating temperature.

DEFINITION: Pixel-to-pixel dark field nonuniformity
This is a repeatable pixel-to-pixel video output level disparity, observed as a quasi-random noise pattern across a frame. It is sourced in individual pixel dark current variations and increases with increasing gain, increasing exposure time and increasing operating temperature. It is sometimes called fixed pattern noise.

DEFINITION: Dark field Vertical Shading
This is a disparity in the mean video output level along the vertical dimension of a frame as compared to mean level for the frame. The disparity takes the form of a slow change in the mean video output level for adjacent horizontal video lines. Dark field Vertical shading may or may not be repeatable from frame-to-frame. It increases with increasing gain, increasing exposure time, and increasing ambient operating temperature.

DEFINITION: Dark field Horizontal shading [DROOP]
This is a disparity in the mean video output level across the horizontal dimension of a frame as compared to mean level for the frame. The disparity takes the form of a slow change in the mean video output level for adjacent vertical video columns. Dark field Horizontal shading is repeatable from frame-to-frame. It increases with increasing gain, increasing exposure time, and increasing ambient operating temperature.

DEFINITION: Dark field Vertical line nonuniformity
This is a disparity in the mean video output level for single or groups of vertical video columns that takes the form of an abrupt change between immediately adjacent or nearby vertical video columns. It is repeatable from frame to frame and increases with increasing gain, increasing exposure time, and increasing/ ambient operating temperature.

DEFINITION: Dark field Horizontal line nonuniformity
This is a disparity in the mean video output level for single or groups of horizontal video lines that takes the form of an abrupt change between immediately adjacent or nearby horizontal video columns. It is repeatable from frame to frame and increases with increasing gain, increasing exposure time, and increasing/ ambient operating temperature.

Model ES-1.0 images exhibit *dark field noise* sourced in the camera sensor and electronics. The noise is substantially random with an approximately Poisson distribution. It is independent of illumination level and exists even in dark field frames captured with little or no illumination. The observed effect is a variation in any pixel's digital video output count from frame-to-frame. Note that this behavior is not illumination-related *shot noise*.

Camera images also exhibit dark field nonuniformities. The nonuniformities are referred to as *pixel-to-pixel dark field nonuniformity*, *vertical line nonuniformity*, *horizontal line nonuniformity*, *vertical shading*, and *droop*. These nonuniformities are independent of illumination level and exist even in dark field frames captured with little or no illumination.

Reported *dark field noise and nonuniformity* measures discount known sensor defects.

FIGURE 4 reports Model ES 1.0 *pixel-to-pixel dark field noise and nonuniformity* for DGN 1-4 digital gain.

FIGURE 5 reports Model ES 1.0 *pixel-to-pixel dark field noise and nonuniformity* for DGN 1-4 digital gain as function of ambient operating temperature.

Pixel-to-pixel dark field noise and nonuniformity is an average for a 10,000 pixel region in the center of the frame, and is reported as percent RMS. of full scale. Full scale is a digital video output count of 255 or 1023 for 8 and 10 bit output respectively.

Dark field horizontal and vertical shading, and horizontal and vertical line nonuniformity measurements are derived from column or row averages over 200 columns or rows in the center of the frame for the entire vertical or horizontal dimension of the frame, and are reported as percent peak of full scale. Full scale is a digital video output count of 255 or 1023 for 8 and 10 bit output respectively.




FIGURE 4




FIGURE 5

Measurement conditions:

*Complete darkness
0.0325 seconds exposure time
Digital video output
Dual channel mode
Defect concealment disabled
ROI image samples*

*Column averaging used for droop
Column averaging used for vertical line
nonuniformity.*

*Row averaging used for vertical shading
Row averaging used for horizontal line
nonuniformity.*

8) Nonlinearity

DEFINITION: Nonlinearity

This is defined as the measure of departure from the camera's theoretical GAMMA TRANSFER FUNCTION. Nonlinearity is computed by applying the camera's inverse theoretical GAMMA TRANSFER FUNCTION to the specified range of X axis coordinate values, and computing the departure of the resultant function from perfect linearity.

GAMMA TRANSFER FUNCTION

The line is defined and fit for the range between the Y axis coordinate value (i.e. mean signal output axis) of the lowest VALID mean signal output level (darkness), and the Y axis coordinate value of the highest VALID mean signal output.

PHOTON TRANSFER FUNCTION

The line is defined and fit for the range between the X axis coordinate value (i.e. mean signal output axis) for the lowest VALID mean signal output level (darkness), and the Y axis coordinate value (i.e. noise variance axis) for highest VALID mean signal output.

A best fit line to the GAMMA TRANSFER FUNCTION is computed as follows: For the function $y = f(x) = ax + b$, where y and x are output gray level and incident illumination respectively on the camera's gamma curve, a and b are calculated so that the sum of the squares of the errors given by $\sum [f(x_i) - y_i]^2$ is minimized. The measure of Nonlinearity may be reported in any or all of three ways.

% RMS. - Expressed as a percent, and defined as the ratio between the average Y axis departure from the theoretical GAMMA TRANSFER FUNCTION and the X axis value for the highest Valid mean signal output.

% INTEGRAL - Expressed as a percent, and defined as the ratio between the normalized absolute value sum of all Y axis departures from the theoretical GAMMA TRANSFER FUNCTION and the X axis value for the highest Valid mean signal output.

% PEAK DIFFERENTIAL - Expressed as a percent, and defined as the ratio between the maximum Y axis departure between any two adjacent X axis coordinate values from the theoretical GAMMA TRANSFER FUNCTION and the X axis coordinate value for the highest Valid signal output.

<u>GAIN</u> (DGN 'X')	<u>GAMMA TRANSFER FUNCTION</u>		<u>PHOTON TRANSFER FUNCTION</u>	
	<u>NONLINEARITY</u> (8 BIT %RMS.)	<u>NONLINEARITY</u> (8 BIT % Integral)	<u>NONLINEARITY</u> (8 BIT %RMS.)	<u>NONLINEARITY</u> (8 BIT % Integral)
DGN 1	0.16	1.5	0.13	0.85
DGN 2	0.24	2.25	0.02	0.02
DGN 4	0.57	5.14	0.04	0.37

<u>GAIN</u> (DGN 'X')	<u>GAMMA TRANSFER FUNCTION</u>		<u>PHOTON TRANSFER FUNCTION</u>	
	<u>NONLINEARITY</u> (10 BIT %RMS.)	<u>NONLINEARITY</u> (10 BIT % Integral)	<u>NONLINEARITY</u> (10 BIT %RMS.)	<u>NONLINEARITY</u> (10 BIT % Integral)
DGN 1	0.48	17.64	0.01	0.12

NOTES: When imaging with very high illumination levels (> 255 or > 1023 for 8 or 10 bit output), the slow onset of sensor blooming may give rise to increased nonlinearity.

When imaging with low high illumination levels and very high gain (< 5 counts and 24 dB), photoelectron traps may also give rise to increased nonlinearity.

Measurement conditions:

Infrared-filtered Xenon light source
 F11 nonimaging illumination

Uncoated clear glass window installed
 0.0325 seconds exposure time
 25 ° C ambient operating temperature
 Digital video output
 Dual channel mode
 Defect concealment disabled
 Center of frame ROI image samples

Labsphere Model ISP 4000 Photometer used for relative illuminance measurements.

9) **Minimum irradiance and illuminance**

DEFINITION: Minimum illumination and Minimum irradiance
 Defined on the GAMMA TRANSFER FUNCTION as the X axis coordinate value (i.e. illuminance or irradiance level) where for the PHOTON TRANSFER FUNCTION, the [SQRT(Y axis coordinate value)] (i.e. signal variance level) and the X axis coordinate value (i.e. signal output level) are equal. In other words, the illumination or irradiance level required to produce a mean video level output equal to the RMS. noise.

Minimum irradiance is reported in nanowatts per cm² @ 550 nm's. Horizontal and vertical shading, and known sensor defects are not considered.

<u>GAIN</u> (DGN 'X')	<u>Minimum irradiance</u> (8 bit)	<u>Minimum irradiance</u> (10 bit)	
1	3.3	2.1	nanowatts per cm ² @550 nm's
2	2.9	NA	"
4	2.7	NA	"

<u>GAIN</u> (DGN 'X')	<u>Minimum illuminance</u> (8 bit)	<u>Minimum illuminance</u> (10 bit)	
1	-	-	millilux for XENON illumination
2	-	-	"
4	-	-	"

10) Photoresponse nonuniformities

The Model ES 1.0 exhibits nonuniformities arising from illumination of the camera's sensor. Photoresponse nonuniformities occur as *pixel-to-pixel photoresponse nonuniformity, vertical shading, droop, vertical line nonuniformity, and horizontal line nonuniformity.*

<p>DEFINITION: Pixel-to-pixel photoresponse nonuniformity Photo-response non-uniformity arises from variations in individual pixel photo-responsivity. The result is differing video level outputs for pixels that are equally illuminated. The observed effect can be individual pixel-to-pixel variations and/or regional variations across a single frame.</p>
<p>DEFINITION: Bright field Vertical shading This is a disparity in the mean video output level across the vertical dimension of a frame as compared to mean level for the frame. The disparity takes the form of a slow change in the mean video output level for adjacent horizontal video lines. Bright field Vertical shading is repeatable from frame-to-frame. Bright field Vertical shading is not substantially influenced by gain, exposure time or temperature.</p>
<p>DEFINITION: Bright field Horizontal shading [DROOP] This is a disparity in the mean video output level across the horizontal dimension of a frame as compared to mean level for the frame. The disparity takes the form of a slow change in the mean video output level for adjacent vertical video columns. Bright field Horizontal shading is uniform across vertical video columns, and is repeatable from frame-to-frame. Bright field Horizontal shading is not substantially influenced by gain, exposure time or temperature.</p>
<p>DEFINITION: Bright field Vertical line nonuniformity This is a disparity in the mean video output level for single or groups of vertical video columns that takes the form of an abrupt change between immediately adjacent or nearby vertical video columns. It is repeatable from frame to frame and not substantially influenced by increasing gain, increasing exposure time, or increasing ambient operating temperature.</p>
<p>DEFINITION: Bright field Horizontal line nonuniformity This is a disparity in the mean video output level for single or groups of horizontal video lines that takes the form of an abrupt change between immediately adjacent or nearby horizontal video columns. It is repeatable from frame to frame and not substantially influenced by increasing gain, increasing exposure time, or increasing ambient operating temperature.</p>

FIGURE 6 reports *pixel-to-pixel photoresponse nonuniformity, vertical shading, droop, and horizontal line nonuniformity* as a percent of the camera's mean gray level output.

Pixel-to-pixel photoresponse nonuniformity is derived from a multiple frame average. It is an average for a 10,000 pixel region in the center of the frame, and reported as percent RMS. of the mean signal output.

Vertical shading, droop, and horizontal line nonuniformity are reported as percent peak.

Bright field vertical shading and droop, and vertical and horizontal line nonuniformity measurements are derived from row or column averages over 200 columns or rows in the center of the frame for the entire vertical or horizontal dimension of the frame, and are reported as percent peak of the mean signal output.

All measures exclude all dark field nonuniformities, exposure nonuniformity, and known sensor defects.




FIGURE 6

NOTES: For the ES 1.0, Droop may be greater at the DGN 1 setting when operating with exposure times of 10 milliseconds or less.

Measurement conditions:

*Infrared-filtered Xenon light source
F11 nonimaging illumination
Uncoated clear glass window installed
0.0325 seconds exposure time
25 ° C ambient operating temperature
Digital video output
Dual channel mode
Defect concealment disabled
Center of frame ROI image samples
Frame averaging used for pixel-to-pixel
photoresponse nonuniformity.*

*Column averaging used for droop
and vertical line nonuniformity .*

*Row averaging used for vertical shading
and horizontal line nonuniformity .*

11) PHOTON TRANSFER FUNCTION and GAMMA TRANSFER FUNCTION

These are the primary MODEL ES 1.0 transfer functions from which many other reported performance measures are derived.

DEFINITION: PHOTON TRANSFER FUNCTION

Defined as a two dimensional discrete function describing the relationship between the *mean* and the *variance of the random component* of the camera's signal output. The X axis is the camera's *mean signal output* and the Y axis is the *variance of the signal output*. The function is defined from complete darkness up through clipping of the A/D where the camera's digital output range is exceeded.

DEFINITION: GAMMA TRANSFER FUNCTION

Defined as a two dimensional discrete function describing the relationship between a camera's *mean signal output* and the *illuminance/irradiance, or luminous/radiant* energy delivered to the camera. The X axis is the *input irradiance* and the Y axis is the camera's *mean signal output*. The plot is from complete darkness up through clipping of the A/D where the camera's digital output range is exceeded.

DEFINITION: ELECTRONIC GAIN

Defined as the *derivative* of the best fit line to the PHOTON TRANSFER FUNCTION. Reported in counts per photoelectron.

NOTES: The QUANTUM EFFICIENCY computed from the PHOTON TRANSFER FUNCTION is for the complete CAMERA SYSTEM, not the sensor alone. It includes other optical losses associated with the sensor cover glass and clear glass window.

FIGURES 7,8,9, and 10 report MODEL ES 1.0 8 and 10 bit PHOTON TRANSFER FUNCTIONS and GAMMA TRANSFER FUNCTIONS for DGN 1,2, and 4 gain settings and exposure time of 0.0325 seconds. Photoresponse and dark field nonuniformity, and black level offset have been subtracted out.

The Model ES 1.0 GAMMA TRANSFER FUNCTIONS are overlaid onto the PHOTON TRANSFER FUNCTIONS. FIGURE 7,8,9, and 10 are four axis plots. The bottom and left axes describe the PHOTON TRANSFER FUNCTIONS. The right and top axes describe the GAMMA TRANSFER FUNCTIONS.

The PHOTON and GAMMA TRANSFER FUNCTIONS for broad band illumination are not reported.

The measurements are an average over a 250,000 pixel region in the center of the frame and exclude dark field nonuniformities, exposure nonuniformity, and known sensor defects.

<u>GAIN</u> (DGN 'X')	<u>ELECTRONIC GAIN</u> (8 bit counts/photoelectron)	<u>GAIN</u> (DGN 'X')	<u>ELECTRONIC GAIN</u> (10 bit counts/photoelectron)
DGN 1	0.006	DGN 1	0.023
DGN 2	0.012	DGN 2	NA
DGN 4	0.026	DGN 4	NA




FIGURE 7




FIGURE 8




FIGURE 9

DGN 1




FIGURE 10

Measurement conditions:

*Infrared-filtered Xenon light source
F11 nonimaging irradiance
Irradiance measured at the lens mount
flange plane.*

*Uncoated clear glass window installed
25 ° C ambient operating temperature
0.0325 seconds exposure time
Digital video output
Dual channel mode
Defect concealment disabled
Center of frame ROI image samples*

*Frame subtraction used to remove
photoresponse nonuniformities.*

12) MODULATION TRANSFER FUNCTION

DEFINITION: MODULATION TRANSFER FUNCTION

The MODULATION TRANSFER FUNCTION (MTF) is the modulation captured by the camera when imaging sine waves of different spatial frequencies. Sine waves of known modulation are imaged into the Camera and modulation is computed by: $(\text{Max.} - \text{Min.}) / (\text{Max.} + \text{Min.})$, as a function of increasing spatial frequency. The MTF may change with the spectrum of the illumination.

MODULATION TRANSFER FUNCTION (MTF) is reported in FIGURE 11. The measurements are in the center of the frame using optimum target phasing. Horizontal and vertical MTF are nearly equal.

A Nikon 95 mm printing lens was measured for MTF and used to image a sine wave test target into the camera. The lens MTF was then divided out to yield the Camera MTF.

Reported measurements are column averages over approximately 25 rows in the center of the frame. Horizontal and vertical shading, all dark field nonuniformities, and known sensor defects are intentionally avoided.




FIGURE 11

NOTES: Average Camera MTF, (i.e. less than optimum target phasing) will yield lower modulation figures for the higher spatial frequencies.

MTF will change with the spectrum of the illumination.

MTF drops rapidly during blooming suppression.

Effects of transfer efficiency, and crosstalk, and linear smear within the camera sensor result in slight variations in MTF depending upon location within the camera frame. The variations will be most noticeable at higher spatial frequencies using: high camera gain settings and low light levels, or low camera gain settings and high light levels.

MTF is measured using 8 bit output.

NOTES: MTF is measured using DGN 1 gain setting.

Measurement conditions:

*Infrared-filtered Xenon light source
 Integrating sphere illuminator

 Uncoated clear glass window installed
 0.0325 seconds exposure time
 DGN 1 gain setting
 25 ° C ambient operating temperature
 Digital video output
 Dual channel mode
 Defect concealment disabled
 Nikon 95 mm printing lens
 F 5.6
 Center of frame ROI image samples
 SINE PATTERNS INC. Sine wave test pattern
 2 times target demagnification*

13) Limiting Resolution

DEFINITION: CONTRAST TRANSFER FUNCTION
 The CONTRAST TRANSFER FUNCTION (MTF) is the contrast captured by the camera when imaging black/white line pairs waves of different spatial frequencies. Black/white line pairs of known contrast are imaged into the Camera and contrast is computed by: $(Max.-Min.)/(Max.+Min.)$, as a function of increasing spatial frequency. The CTF may change with the spectrum of the illumination. The CTF of the imaging lens is compounded into the measure and NOT subtracted out.

DEFINITION: Limiting resolution
 The Limiting resolution is the maximum imaged spatial frequency, reported in cycles or TV lines/picture height, wherein the camera can maintain a minimum contrast value of 30%.

<u>Illumination</u>	<u>Gain</u> (DGN 'X')	<u>Limiting Resolution</u> (cycles/picture height)
Infrared-filtered Xenon	1	500
	(dB)	(TV lines/picture height)
Infrared-filtered Xenon	1	-

NOTES: 1 cycle is equivalent to 2 TV lines using the SMPTE standard for measurement of resolution. Also, the measured contrast for TV lines will be different (usually higher) than the measured modulation using sine waves.

14) CONTRAST TRANSFER FUNCTION

DEFINITION: CONTRAST TRANSFER FUNCTION

The CONTRAST TRANSFER FUNCTION (CTF) is the contrast captured by the camera when imaging black/white line pairs waves of different spatial frequencies. Black/white line pairs of known contrast are imaged into the Camera and contrast is computed by: $(\text{Max.} - \text{Min.}) / (\text{Max.} + \text{Min.})$, as a function of increasing spatial frequency. The CTF may change with the spectrum of the illumination. The CTF of the imaging lens is compounded into the measure and **NOT** subtracted out.

The Camera CONTRAST TRANSFER FUNCTION (CTF) is reported in FIGURE 12. The measurements are in the center of the frame using optimum target phasing. Horizontal and vertical CTF are nearly equal.

A Nikon 60 mm AF lens is used to image the line pair test target into the camera. The CTF of the Nikon lens is compounded into the measure and **NOT** divided out.

Reported measurements are column averages over approximately 25 rows in the center of the frame. Horizontal and vertical shading, all dark field nonuniformities, and known sensor defects are intentionally avoided.

NOTES:

Average Camera CTF, (i.e. less than optimum target phasing) will yield lower contrast figures for the higher spatial frequencies.

CTF will change with the spectrum of the illumination.

CTF drops rapidly during blooming suppression.

Effects of transfer efficiency, and crosstalk within the camera sensor result in slight variations in CTF depending upon location within the camera frame. The variations will be most noticeable at higher spatial frequencies using high camera gain settings and low light levels or low camera gain settings and high light levels.

CTF is measured using 8 bit output.




FIGURE 12

Measurement conditions:

Infrared-filtered Xenon light source

*Uncoated clear glass window installed
 0.0325 seconds exposure time*

25 ° C ambient operating temperature

Digital video output

Dual channel mode

Defect concealment disabled

Nikon 60 mm AF lens

F 5.6

Center of frame ROI image samples

NBS type microscopy test target

1 times target demagnification

15) Spectral responsivity

DEFINITION: Spectral responsivity
 Spectral responsivity is defined as the Camera's responsivity when imaging with different spectral wavelength bands.

FIGURE 13 reports Model ES 1.0 spectral responsivity in counts per uw/cm²-second rounded to the nearest 10 counts. Reported responsivity is for the camera exclusive of any filters, and in conjunction with the Schott infrared filter option. When used with lenses or other filters the net spectral response may be different.




FIGURE 13

SPECTRAL RESPONSIVITY

(counts/uw/cm²/second)

<u>WAVELENGTH</u>	<u>DGN 1 [8 bits]</u>	<u>DGN 2 [8 bits]</u>	<u>DGN 4 [8 bits]</u>	<u>DGN 1 [10 bits]</u>
400	3550	6820	13460	14780
450	3880	7460	14710	16150
500	4690	9020	17780	19530
550	5020	9650	19030	20900
600	4860	9360	18450	20270
650	4450	8550	16870	18530
700	2660	5120	10090	11080
750	1900	3660	7210	7920
800	1420	2730	5380	5910

<u>WAVELENGTH</u>	<u>DGN 1 [8 bits] SCHOTT</u>	<u>DGN 2 [8 bits] SCHOTT</u>	<u>DGN 4 [8 bits] SCHOTT</u>	<u>DGN 1 [10 bits] SCHOTT</u>
400	3130	6030	11880	13050
450	3400	6540	12900	14170
500	4190	8070	15900	17470
550	4490	8630	17020	18690
600	4210	8100	15980	17550
650	3390	6510	12850	14110
700	11570	3030	5970	6560
750	760	1460	2880	3170
800	350	670	1320	1450

Measurement method:

*Monochromator (per Kodak MTD sensor data)
 Measured at the lens mount flange plane*

Measurement conditions:

*Uncoated clear glass window installed
 25 ° C ambient operating temperature
 0.0325 seconds exposure time*

Cross-check at 550 nm:

EG&G Model DR 2550 Radiometer

16) **Crosstalk and Smear**

DEFINITION: Crosstalk

This is defined as either or both of the following two effects:

- 1) Photocurrent integration in pixels that are not illuminated, but are nearby to illuminated pixels.
- 2) Photocurrent integration in pixel transfer registers resulting from light leakage through the sensor substrate. This results in premature blooming at the pixels of entry, along with nonlinear GAMMA TRANSFER FUNCTION behavior and premature blooming at other pixel locations.

Crosstalk increases linearly with increasing exposure time, higher illumination level, and increasing illumination wavelength.

DEFINITION: Smear

This is defined as an apparent linear smear of an imaged object caused by photocurrent integration in transfer registers (effect #2). Smear increases linearly with the luminous exposure time until the entire vertical dimension of the frame is affected.

FIGURES 14 and 15 reports the effect of MODEL ES 1.0 *crosstalk* arising from photocurrent integration in vertical pixel transfer registers. *Crosstalk* is reported as percent of the mean digital video output count. Reported measures are for the DGN 2 gain setting and the worst case pixel; namely, that pixel whose video output level receives the greatest *crosstalk* contribution from other vertically adjacent pixels.

FIGURE 16 reports MODEL ES 1.0 *smear* resulting from *crosstalk*.




FIGURE 14

NOTES: The crosstalk may be different for different DGN gain settings.

Measurement conditions:

*DGN 2 gain setting
 Infrared-filtered Xenon light source
 Continuous wave illumination*

*Uncoated clear glass window installed
 25 ° C ambient operating temperature
 Center of frame ROI image samples
 Digital video output
 Dual channel mode
 Defect concealment disabled*




FIGURE 15

NOTES: The crosstalk may be different for different DGN gain settings.

Measurement conditions:

*DGN 2 gain setting
 Infrared-filtered Xenon light source
 Continuous wave illumination*

*Uncoated clear glass window installed
 25 ° C ambient operating temperature
 Center of frame ROI image samples
 Digital video output
 Dual channel mode
 Defect concealment disabled*




FIGURE 14

Measurement conditions:

Infrared-filtered Xenon light source

*Uncoated clear glass window installed
Back illuminated variable height
horizontal slit.*

*25 ° C ambient operating temperature
Center of frame ROI image samples
Digital video output
Dual channel mode
Defect concealment disabled*

17) Blooming suppression

DEFINITION: Over-illumination

This is defined as an imager illumination level that is higher than the maximum value of the camera's digital count output range. Illumination levels higher than this maximum digital count cannot be 'sensed' by the Camera.

DEFINITION: Blooming

This is defined as a dynamic 'pixel overflow' condition within many sensors, usually caused by over-illumination. Pixels are illuminated beyond their ability to integrate charge, causing charge overflow into adjacent structures. Adjacent structures may be nearby pixels, transfer registers, or other sensor structures. The observed effect is usually a 'growing vertical stripe or bright spot', or other bright and/or distorted/extended image objects that appear to have been imaged abnormally by the camera. The resolution of an imager is reduced rapidly and significantly as blooming begins to occur, and can reach zero depending on the specific region of the camera sensor. An important identifying characteristic of blooming is its rapid spreading across the entire vertical dimension of the frame.

DEFINITION: Blooming suppression factor

This is defined as the camera over-illumination multiplication factor (x) wherein the Camera *will not bloom*.

Imager *blooming* is observed as a *growing bright spot* which expands more rapidly in the vertical axis than the horizontal, eventually resulting in a vertical stripe. *Over-illumination* of the camera may produce this condition.

FIGURE 17 reports MODEL ES 1.0 *blooming suppression factor* for the DGN 2 gain setting.




FIGURE 15

NOTES: The blooming suppression factor may be different for different DGN gain settings.

Measurement conditions:

*DGN 2 gain setting
 Infrared-filtered Xenon light source
 Continuous wave illumination*

*Uncoated clear glass window installed
 25 ° C ambient temperature
 Center of frame ROI image samples
 Digital video output
 Dual channel mode
 Defect concealment disabled*

18) Stray Light

19) Valid signal output and measurement accuracy

DEFINITION: Valid signal output
 Valid signal output is defined as that portion of the GAMMA and PHOTON TRANSFER FUNCTIONS where the reported performance figure is determined to be a '99% confident' measure of the camera's performance.
 'Confidence' is defined as the standard statistical 99% confidence measure of a parameter of a population.
 The parameters for which the confidence applies are the reported MEAN and VARIANCE figures of the GAMMA and PHOTON TRANSFER FUNCTIONS.
 The population for which the confidence applies is the population of measures or estimates used in reporting the MEAN and VARIANCE figures of the GAMMA and PHOTON TRANSFER FUNCTIONS.

<u>MEASURE</u>		<u>ACCURACY</u>	<u>UNITS</u>
PHOTON TRANSFER FUNCTION		+/- 1	percent of value
GAMMA TRANSFER FUNCTION		+/- 1	percent of value
MODULATION TRANSFER FUNCTION	Relative	+/- 5	percent of value
	Absolute	+/- 15	"
CONTRAST TRANSFER FUNCTION	Relative	+/- 5	percent of value
	Absolute	+/- 15	"
Responsivity		+/- 15	percent of value
Gamma ratio		+/- 1	percent of value
Blooming suppression		+/- 10	percent of value
Crosstalk		+/- 10	percent of value
Smear		+/- 10	percent of value
Nonlinearity		+/- 1	percent of value
Dynamic range		+/- 1	percent of value
Dark Field noise and nonuniformity		+/- 1	percent of value
Maximum Luminous signal/noise ratio		+/- 1	percent of value
Minimum illuminance and irradiance		+/- 10	percent of value
Photoresponse nonuniformity		+/- 1	percent of value
Stray light		-	percent of value
Limiting Resolution	Relative	+/- 5	percent of value
	Absolute	+/- 15	"
Exposure variability and nonuniformity	Relative	+/- 1	percent of value
	Absolute	+/- 1	"
Spectral responsivity	Consult Eastman Kodak Company Motion Analysis Systems Division		
Equivalent ISO speed	Computed from photopic responsivity measure		

20) Measurement Instrumentation and data reduction tools

INSTRUMENTATION

ILC 300 watt xenon lamp and linear power supply
Oriental cold mirror and liquid light guide
Labsphere 6" integrating sphere with 2 " port
Labsphere 4" integrating sphere with 1 " port
Kodak custom F11 baffle tube
Sine Patterns sine wave target
Melles Griot 10nm band width 550 nm filter
Labsphere model ISP4000 integrating sphere photometer
EG&G Model 2550 Radiometer/photometer
Nikon 95 MM Printing Lens
Nikon 60 mm photographic lens
Precision Digital Images Model AD4 digital frame grabber
Macintosh IIfx computer system

DATA REDUCTION TOOLS

"Image 1.55.3_VDM" image processing software
"Image 1.58_VDM" image processing software
"Deltagraph" data analysis and plotting software
"Statview" statistical data analysis software

NOTES

Redlake / Kodak Megaplug ES 310 Summary

MAX: PCI-1422 img0 684 x 242 Disable serial

Note the use of modified camera file "Redlake ES 310 Modified"
in NI-IMAQData.

Serial communications via COM1

TTermPro: 9600 baud / 8 data / 1 stop / no parity
Configure TTermPro->Terminal for local echo

Issue the following serial commands (followed by ENTER)

IDN? (camera responds with ID string)
VID OF (disable analog video output, unless needed)
BLK ON (block mode)
BST 1 (block start line 1)
BSP 121 (required for 684 W x 242 H frame)
FRS 60 (set camera frame rate to 60 FPS [85 max])
AEX ON (autoexposure mode)

Run the LabVIEW application TCP VIDEO

nframes specifies the size of the circular buffer (200 frames)
initialize "trialname" to an appropriate string identifier

The application will create a folder on the desktop with the
name entered in the trialname box. Binary dumps of the circular
buffer are named sequentially: 000.vid, 001.vid, . . . and saved
to this folder.

Video files can be replayed with the MATLAB application "playback".

From another PC:

```
telnet  
set crlf ( terminator sequence CR/LF appended to messages )  
set localecho ( TCP VIDEO does not echo )  
set escape x ( if necessary )  
open 132.239.185.132 6342
```

```
send "s" to write the circular buffer to disk  
send "q" to cause TCP VIDEO to exit
```

Notes

The cable connections are fragile and often result in line dropout.
See additional camera documentation in lab files.