

HISTOSLIDE 2000

A compact sliding microtome







Precision engineering

The newly designed and virtually maintenance-free micrometer feed system is totally enclosed to prevent section waste or cleaning fluids from entering the mechanism, where accumulated debris will eventually cause failure.

Smooth, resistant surfaces make the HISTOSLIDE 2000 easy to clean. For additional convenience, an optional waste tray is available.




The section selector and the coarse feed are positioned conveniently on the side of the instrument within easy reach of the user. The selected section thickness can be read from the scale at a glance.




Increased productivity

The specimen advance can be either automatically or manually initiated. When sectioning smaller blocks, the automatic stop can be easily repositioned thus shortening the stroke saving both time and effort. No special tools are required.



Tireless operation

The sledge glides smoothly on five bearings which have been carefully dimensioned to redistribute oil on the slideways. This will minimize the cutting effort without compromising the stability, which is essential for good quality sectioning.




Easy to handle knife holder

The knife can be clamped quickly and is securely held by two screws. The declination can be adjusted by simply loosening the handle — no tools are required. The clearance angle can be adjusted without loosening the retaining screws on the knife. Safety is essential and the knife guard will protect a 22 cm knife.


Rugged specimen holder

The specimen holder will accommodate blocks up to 80×60 mm. Orientation is standard. No tools are required to switch to cassettes. A simple spring loaded



insert is placed in the vise, and sectioning can start immediately.

Reichert-Jung offer the widest range of products for microtomy. All critical components are developed and manufactured at the Reichert-Jung factory in Nussloch



near Heidelberg and are subjected to the strictest quality control to ensure they conform to the highest standards. It is this attention to detail that has established Reichert-Jung as the world's leading manufacturer of microtomes.

Specification and ordering information

Basic specification	Catalog Number
Totally enclosed micrometer system with orientating fixed jaw specimen	
clamp for paraffin blocks. Maximum specimen size 80 × 60 mm. Section	
thickness adjustable in the range of $0-60~\mu m$. Automatic advance in the	
range of $0-30~\mu m$. Selectable in $1~\mu m$ increments from $0-10~\mu m$,	
in 2μ m increments from $10-20 \mu$ m, in 5μ m increments from $20-60 \mu$ m.	
Feed range 40 mm.	
Knife holder for standard knives with knife guard, inclination and	
declination adjustable without loosening the knife.	
Standard system including protective cover	0413_18800
System including oil for slideways 404, micrometer 405 and rust preventative oil 406	90720
Microtome accessories	
Oil for slideways 404, 50 ml	0336 06080
Oil for micrometer 405, 50 ml	0336 06086
Rust preventative oil 406, 100 ml	0336 06093
Waste tray	0413 18960
Cassette clamp insert	0413 18963
Distance plate for knife holder for use with Reichert-Jung cassette clamp 0394 09151	
and freezing equipment such as Frigomobil	0413 18967
Sectioning accessories	
Knife 16 cm profile c (both sides plane)	0216.07100
Knife 16 cm profile d (both sides plane with steep facet)	0216 07132
Knife 22 cm profile c (both sides plane)	0216 07116
Variable knife case for 1 or 2 knives (10, 12 or 16 cm)	0213 11140
Variable knife case for 1 or 2 knives (22, 25 or 30 cm)	0213 11141
Dispenser with 50 disposable blades 75×14 mm (high profile)	0358 12881
Dispenser with 50 disposable blades 75×8 mm (low profile)	0358 13583
Holder for disposable blades (high profile and low profile)	0368 09426
Equipment for frozen sections	
Frigomobil O for specimen cooling	0354 09723
Frigomobil OM for specimen and knife cooling	0354 09717
CO ₂ freezing stage with accessories	0274 05918
Technical data	
Disconsisses	

Dimensions

Width: 30 cm, depth 44 cm, height 33 cm

Weight: 32 kg

Due to a policy of continuous development, we reserve the right to change specifications without notice.



CAMBRIDGE-INSTRUMENTS is an international high-technology group specializing in the development and manufacture of scientific instrumentation, electron optical microscopes, vision care apparatus and semi-conductor manufacturing equipment. With origins that go back over 100 years, Cambridge has an excellent reputation for technical innovation and has pioneered numerous commercial applications of scientific discoveries. Research and development is carried out in the United States and Europe — in Vienna, Heidelberg and Cambridge. In Heidelberg and Vienna, the Reichert-Jung products are produced.

REICHERT-JUNG is the name associated with the world's leading manufacturer of microtomes, ultramicrotomes and research microscopes, which are used extensively in life sciences, histology, clinical pathology, materials science and industry. Over 100 dedicated years of precision engineering experience culminate in the widest range of products for microtomy. The continual investment in research, development and manufacturing facilities are proof of Reichert-Jung's commitment to quality. Four generations of biologists, scientists and technicians have used and are continuing to use Reichert-Jung microtomes both for research and routine applications.

M. O. C. INC.
MICROSCOPICAL OPTICAL CONSULTING
P. O. BOX 586
VALLEY COTTAGE, N.Y. 10989
914-268-6450

Cambridge Instruments Inc.

Sales: P.O. Box 123, Buffalo, New York 14240, U.S.A., Telephone: (716) 891-3000, Tel.: 91-285, Telef. (716) 891-3155 Service: 111 Deer Lake Road, Deerfield, Illinois 60015, U.S.A., Telephone: (312) 405-0123, Telefax: (312) 405-0147

MOC

Microscopical
Optical
Consulting inc.

P.O Box 586 Valley Cottage, New York 10989 (914) 268-6450 FAX 914-268-0172

DATE FCO. 1, 1994

TO: DAVID KIEINSELD

FAX: 908/582-2451

FROM: DEALLA R. FAGANO

PHONE NO. 914/268-6450

FAX: 914/268-0172

2 PAGES INCLUDING COVER SHEET.

comments: I added to the gutte unat is reloted & what you wanted.

I hope this is okay.

SINCERELY YOURS, questions please calla.

M.O.C. INC.

Dearra L. Lagaro

MOC

Microscopical Optical Consulting inc.

P.O. Box 586, Valley Cottage, NY 10989 (914) 268-6450 / FAX 914-268-0172 QUOTATION

ORDER NUMBER

Nº 7054

02-01-94

Mr. David Kleinseld ROOM 1C-463 AT&T BELL LABORATORIES 600 Mountain Avenue Murray Hill, NJ 07974

Fax#: 908/582-2451

PRICES AND TECHNICAL SPECIFICATIONS HEREIN ARE SUBJECT TO CHANGE WITHOUT NOTICE.

TERMS		DELIVERY	SHIP	VIA	
NET 30 DAYS					
CATALOG NUMBER			. DESCRIPTION	•	
045325121	1	mection thickness Maintenance free rollers with tota system and orient clamp for paraffi for standard kniv blade holders. I adjustable withou installed guide rentering the rolle advance in the rai increments from 0- 10-20um, in 5um in	complete with Autoradjustment with: Sliding microtome of the control of the contr	n guide eter feed cimen holder disposable instion angles Vertically is from atic specimen ctable in lum ments from Oum. Total tive cover	6,936.00
027405918 021404121 021607100 021607132	1 1 1 1	Cooling Stage Bottle Stand 16cm C Profile Kn 16cm D Profile Kn SHIPPING CHARGES 35-45 DAYS DELIVER YOUR PURCHASE ORDE	ife \$ 140.00 RY TIME, AFTER RECE	\$ \$ \$ \$ IPT OF	879.00 637.00 216.00 216.00
·		PRICE QUOTATION VA	ALID FOR (30) DAYS.		

This equipment is subject to modifications and the above quotation does not constitute a commitment to sell these items.

All prices are subject to local or state sales taxes, if any.

Authorized Signature Least & Mark

MOC

Microscopical Optical Consulting Inc.

P.O. Box 586, Valley Cottage, NY 10989 (914) 268-6450 / FAX 914-268-0172 QUOTATION

ORDER NUMBER 7054

02-01-94

Mr. David Kleinseld ROOM 1C-463 AT&T BRLL LABORATORIES 600 Mountain Avenue Murray Hill, NJ 07974

Fax#: 908/582-2451

PRICES AND TECHNICAL SPECIFICATIONS HEREIN ARE SUBJECT TO CHANGE WITHOUT NOTICE.

TERMS		DELIVERY		SHIP VIA	·-·
NET 30 DAYS			:		
CATALOG NUMBER			DESCRIPTION		
			Y		
HISTOSLIDE 2000 R complete with Automatic section thickness adjustment with: Maintenance free Sliding microtome on guide rollers with totally enclosed micrometer feed system and orientating fixed-jaw specimen clamp for paraffin sectioning. Knife holder for standard knives and all current disposable blade holders. Inclination and declination angles adjustable without loosening knife. Vertically installed guide rollers prevent debris from entering the roller bearings. Automatic specimen advance in the range of 1-30um, selectable in lumincrements from 0-10um, in 2um increments from 10-20um, in 5um increments from 20-60um. Total specimen feed range is 40mm. Protective cover included. Factodes puttur call stage.		6,936.00			
021607100		16cm Knife 1821	slack)		
		SHIPPING CHARGES	\$ 140.00		
		35-45 DAYS DELIVE YOUR PURCHASE ORD		RECEIPT OF	
		PRICE QUOTATION V	ALID FOR (30) DA	YS.	
				• .	
	- 1 1		•	•	

This equipment is subject to modifications and the above quotation does not constitute a commitment to self these liems.

All prices are subject to local or state sales taxes, if any.

M Microscopical O Optical C Consulting inc.

P.O. Box 586, Valley Cottage, NY 10989 (914) 268-6450 / FAX 914-268-0172

QUOTATION

ORDER NUMBER

Nº 7054

DATE

02-01-94

Mr. David Kleinseld ROOM 1C-463 AT&T BELL LABORATORIES 600 Mountain Avenue Murray Hill, NJ 07974

Fax#: 908/582-2451

PRICES AND TECHNICAL SPECIFICATIONS HEREIN ARE SUBJECT TO CHANGE WITHOUT NOTICE.

TERMS		DELIVERY			SHIP VIA		
NET 30 DAYS			•				
CATALOG NUMBER				DESCRIPTION	· · · · · · · · · · · · · · · · · · ·	·	·
045325121	1	section the Maintenance rollers wire system and clamp for	th totally e orientating paraffin sec	stment wi ng microt enclosed m fixed-ja tioning.	th: ome on guide icrometer fe		
		adjustable installed entering t advance in increments 10-20um, i	without loc guide roller he roller be the range of from 0-10um	sening knes prevent arings. of 1-30um, n, in 2um aents from	declination ife. Vertice debris from Automatic spectable increments f 20-60um. To	ally ecimen in lum rom tal ver	3,936.00
027405918 021404121 021607100 021607132	1 1 1 1	16cm D Pro	and ofile Knife ofile Knife CHARGES \$ 14	40.00 ME, AFTER	RECEIPT OF	\$ \$ \$ \$	879.00 637.00 216.00 216.00
		PRICE QUOT	ATION VALID	FOR (30)	DAYS.		

This equipment is subject to modifications and the above quotation does not constitute a commitment to sell these items.

All prices are subject to local or state sales taxes, if any.

Authorized Signature

Canak Hogano

Microscopical Optical Consulting inc.

P.O. Box 586, Valley Cottage, NY 10989 (914) 268-6450 / FAX 914-268-0172

Q	UO	TAT	1OI

ORDER NUMBER
102 7054

DATE 02-01-91

Mr. David Kleinseld ROOM 1C-463 AT&T BELL LABORATORIES 600 Mountain Avenue Murray Hill, NJ:07974

Fax#: 908/582-2451

PRICES AND TECHNICAL SPECIFICATIONS HEREIN ARE SUBJECT TO CHANGE WITHOUT NOTICE.

TERMS		DELIVERY	SHIP VIA
NET 30 DAYS			·
CATALOG NUMBER	4, .	DESCRIPTION CONTROL OF THE DESCRIPTION CONTROL O	And the second of the second o
045325121	1	HISTOSLIDE 2000 R complete with a section thickness adjustment with Maintenance free Sliding microton rollers with totally enclosed mid system and orientating fixed-jaw clamp for paraffin sectioning. Kr for standard knives and all curre blade holders. Inclination and adjustable without loosening knift installed guide rollers prevent entering the roller bearings. At advance in the range of 1-30um, sincrements from 0-10um, in 2um in 10-20um, in 5um increments from 2 specimen feed range is 40mm. Proincluded.	n: ne on guide crometer feed specimen nife holder ent disposable declination angles fe. Vertically debris from utomatic specimen selectable in lum ncrements from 20-60um. Total
027405918 021404121 021607100 021607132	1 1 1 1	Cooling Stage Bottle Stand 16cm C Profile Knife 16cm D Profile Knife SHIPPING CHARGES \$ 140.00 35-45 DAYS DELIVERY TIME, AFTER F YOUR PURCHASE ORDER. PRICE QUOTATION VALID FOR (30) DA	

This equipment is subject to modifications and the above quotation does not constitute a commitment to sell these items.

All prices are subject to local or state sales taxes, if any.

Authorized Signature

arra K. Fayano

M Microscopical O Optical C Consulting inc.

Sales & Services
Microscopes • Microtomes • Cameras

Helmut Patzig P.O. Box 586 Valley Cottage, NY 10989

(914) 268-6450 FAX 914-268-0172