



Filter set 47 HE
489047 - 9901

excitation: BP 436/25
 beamsplitter: FT 455
 emission: BP 480/40



Filter set 38 HE
489038 - 9901

excitation: BP 470 / 40
 beamsplitter: FT 495
 emission: BP 525 / 50




Filter set 64 HE

489064 - 0000

excitation: BP 587/25
beamsplitter: FT 605
emission: BP 647/70



(typical curves)



Filter set 50

488050 - 9901

excitation: BP 640/30
beamsplitter: FT 660
emission: BP 690/50



(typical curves)

